

Auf der Gröb 2 c
83064 Raubling
Germany
Tel.: +49 (0)8035-9580448
Fax: +49 (0)8035-9580446
info@hernik.de • www.hernik.de

invoice / delivery address

Kom.:

Order

Offer

Helmet CAPOVARIO

Exportprices **Release: 2023**

BASIC

technical specification:

- solid, anatomic frame made of synthetical, size adjustable
- molded rigid foam in 3D-concavity
- microfiber: washable by 40 °C, can be disinfected, breathable, flame resistant, antistatic, dirt- and water-repellent
- forehead protection included
- upholstery inlay (microfiber),
- weight 280 - 320 g
- all accessories can be adapted later, too!



CAPOVARIO

(60050Y-Z), [€ 695,-]



CAPOVARIO LIGHT

(60055Y-Z), [€ 675,-]

size (Y)	0	1	2	3	4
circumference of the head	44 - 48 cm	48 - 51 cm	52 - 56 cm	56 - 60 cm	58 - 62 cm
weight	270 g	280 g	290 g	305 g	320 g
PINKRED	<input type="checkbox"/> 6005X0-P	<input type="checkbox"/> 6005X1-P			
TURQUOISE BLUE	<input type="checkbox"/> 6005X0-T	<input type="checkbox"/> 6005X1-T			
RED			<input type="checkbox"/> 6005X2-R	<input type="checkbox"/> 6005X3-R	<input type="checkbox"/> 6005X4-R
BLACK			<input type="checkbox"/> 6005X2-S	<input type="checkbox"/> 6005X3-S	<input type="checkbox"/> 6005X4-S
GREY			<input type="checkbox"/> 6005X2-G	<input type="checkbox"/> 6005X3-G	<input type="checkbox"/> 6005X4-G
BLUE			<input type="checkbox"/> 6005X2-B	<input type="checkbox"/> 6005X3-B	<input type="checkbox"/> 6005X4-B
strap <small>Please choose ONE of them!</small>	<input type="checkbox"/> standard (600530) [€ 65,-]		<input type="checkbox"/> padded, for size 0, 1, 2 (600535-1) <input type="checkbox"/> padded, for size 3, 4 (600535-2) [€ 115,-]		

X=0 (CV Standard); X=5 (CV LIGHT)


colours:



PINKRED



RED



BLACK



GREY



BLUE




TURQ. BLUE

accessories

4 straps
(600590-X)
[€ 225,-]


size 0
 size 1
 size 2
 size 3
 size 4




closed top side

soft (600520-X) [€ 275,-] hard (600525-X) [€ 545,-]

size 0 size 0
 size 1 size 1
 size 2 size 2
 size 3 size 3
 size 4 size 4




FIDLOCK One-Hand-Use
(600585-X)
[€ 95,-]




upholstery for the buckle
(600580-X)
[€ 70,-]


1: for standard strap
 2: for padded strap



protects the skin - no contact with the buckle

magnetic lock
(600550)
[€ 115,-]


1: for standard strap
 2: for padded strap



note:
• user is not able to put away the helmet
• a special magnetic key must be used!

chin protector
(600515-X)
[€ 155,-]

1: for size 0/1/2
 2: for size 3/4



ear protectors*
(600510-X)
[€ 155,-]

1: for size 0/1/2
 2: for size 3/4



*padded strap necessary!

visor #
(600570)
[€ 475,-]

1: for size 0/1/2
 2: for size 3/4



delivery time 2 - 3 weeks

cheek protector
(600565-X)
[€ 325,-]


1: for size 0/1/2
 2: for size 3/4




spare parts


additional strap, padded
(600535-X)
[€ 115,-]

1: for size 0/1/2
 2: for size 3/4



alternating inlay, 2 straps
(600540-X)
[€ 265,-]

size 0
 size 1
 size 2
 size 3
 size 4



alternating inlay, 4 straps
(600595-X)
[€ 325,-]

size 0
 size 1
 size 2
 size 3
 size 4



production in two steps
(600545-X)
[€ 185,-]

NEW!

A blanks can be adjusted with hot air by a technician to adapt it on place to the individual requirements of the patient's individual requirements before finally obtained and finished in our workshop.

Fixation of all movable elements
(600555-X)
[€ 255,-]


NEW!

After a successful fitting on place, all removable parts are fixed in our workshop. This may be important for very vigorous patients.

Washing then only by hand!

mouth protection, adjustable*
(600560-X)
[€ 395,-]

1: for size 1/2
 2: for size 3/4



*padded strap necessary!