

invoice / delivery address

Reference:

**Order**

**Offer**

# Rehab Car Seat *IPAI* **LGT**

**Exportprices** **Release: 2023**

**BASIC**

**common information (standard):**

- crash tested according to ECE-R44/04
- 5-point-harness with upholstery, height adjustable, can be tightened independently
- large pad between the buckle and the pelvis
- ISOFIX adapter included with one-hand-adjustment
- head rest height adjustable, including belt guides
- Side impact protection SICT, separately adjustable per side
- side supports, individually height-adjustable
- removable, breathable cover (washable at 30 ° C)
- every accessory can be mounted in seconds, even later!



**Details:**

IPAI-LGT	
SEAT WIDTH (pelvis)	31 cm
SEAT WIDTH (knee)	34 cm
SEAT DEPTH	30 cm
BACK HEIGHT	41 - 54 cm
CROWN HEIGHT	59 - 74 cm
HEIGHT LATERALS	20 - 33 cm
DISTANCE LATERALS	14 - 28 cm
angle	85° - 160°
total height	68 - 88 cm
total width	51 cm
total depth	42 cm
weight	ab 7,0 kg
capacity	36 kg

**IPAI-LGT black**  
(460000-E)  
**€ 1.490,00**

**fittings for stabilization**

abduction block, removable

(460400-E)  
[€ 330,-]



head fixation

(460630-E)  
[€ 265,-]



security tray

(460600-E)  
[€ 465,-]



**further options**

protection mat

(to protect the original seats of the car)

(460780-E)  
[€ 80,-]



belt to connect the chest straps

(to improve the position of the belt in front of the breast)

(460640-E)  
[€ 115,-]



angular adjustment (10°)

(460100-E)  
[€ 185,-]



seatdepth adapter

(460850-E)  
[€ 285,-]



**notes**

.....  
.....  
.....  
.....  
.....  
.....  
.....  
.....  
.....