

invoice / delivery address

Reference:

Order

Offer

Rehab Car Seat *IPAI* ^{LGTT}

Exportprices **Release: 2023**

BASIC

common information (standard):

- crash tested according to ECE-R44/04
- 5-point-harness with upholstery, height adjustable, can be tightened independently
- large pad between the buckle and the pelvis
- ISOFIX adapter included with one-hand-adjustment
- head rest height adjustable, including belt guides
- Side impact protection SICT, separately adjustable per side
- side supports, individually height-adjustable
- removable, breathable cover (washable at 30 ° C)
- every accessory can be mounted in seconds, even later!



Details:

IPAI-LGT	
SEAT WIDTH (pelvis)	31 cm
SEAT WIDTH (knee)	34 cm
SEAT DEPTH	30 cm
BACK HEIGHT	41 - 54 cm
CROWN HEIGHT	59 - 74 cm
HEIGHT LATERALS	20 - 33 cm
DISTANCE LATERALS	14 - 28 cm
angle	85° - 160°
total height	68 - 88 cm
total width	51 cm
total depth	42 cm
weight	ab 7,0 kg
capacity	36 kg

**IPAI-LGT
black**
(460000-E)

