

invoice / delivery address

Reference:

☐ **Order**
☐ **Offer**

Rehab Car Seat **STARLIGHT**

Exportprices

Release: 2023

BASIC

Details:

| | STARLIGHT-NXT (standard) | STARLIGHT-NXT (XL-Kit) |
|---------------------|------------------------------|---------------------------|
| SEAT WIDTH (pelvis) | 30 cm | 30 cm |
| SEAT WIDTH (knee) | 35 cm | 35 cm |
| SEAT DEPTH | 26 cm | 32 cm |
| BACK HEIGHT | 27 - 39 cm | 44 - 53 cm |
| CROWN HEIGHT | 51 - 62 cm | 67 - 76 cm |
| HEIGHT LATERALS | 17 - 27 cm | 25 - 38 cm |
| DISTANCE LATERALS | 14 - 19 cm | |
| angle | 90° | |
| total height | 60 - 85 cm | |
| total width | 50 | |
| total depth | 53 | |
| weight | 8,5 kg (without accessories) | |
| capacity | 9 - 36 kg | |

common information (standard):

- crash tested according to ECE-R44/04
- 5-point-harness with upholstery, height adjustable
- belt mechanics installed into the front of the seat
- head rest height adjustable, including belt guides
- additional upholstery (20 mm/40 mm) for the back & headrest
- fold-out cup holder



ATTENTION: ISOFIX adapter or clamp bow necessary for children with poor trunk control or spasticity. Tilting danger!!!



STARLIGHT-NXT
red-black
(480000-E1)



STARLIGHT-NXT
grey-black
(480000-E2)



accessory that makes the base board **NECESSARY**

Necessary for the options below!
(480000-BE)



BASE BOARD



swivel base

☐ Gr. 1 - 3
(480200-E)

☐ right

☐ left



fastener

(independent of the direction of rotation)

☐ inside
(standard)

☐ outside
(480220-E)


☐ angular adjustment (0° - 10°)
(480100-E)

☐ extra angular adjustment + 10°
(10° - 20°)
(480110-E)



☐ ISOFIX adapter
(480575-E)

☐ clamp bow
(480550-E)



☐ Sandwich-Base (mounted under the base board)
(480790-E)



Useful in the following situations:

- with very sporty shaped car seats, in order to compensate for the recess there. This results in a more stable position and more freedom of movement when using the swivel adapter.
- with car Isofix brackets that are very far inside to allow nevertheless contact with the car's backrest and the IPAI-NXT (the ISOFIX adapter sits in the sandwich base)
- approx. 1.5 kg additional weight

foot rest

☐ (480450-E)




☐ swing up*
(480452-E)



*480450-E necessary!


☐ outrigger*
(480453-E)



*480450-E necessary!

foot roll

☐ (480455-E)



fittings for stabilization

head fixation

☐ (480630-E)



security tray

☐ (480600-E)


longer 5-point-harness (+ 15 cm)

☐ (480320-E)


side supports, pair, individually adjustable in height, deep and width

☐ (480800-E)



lateral support for the seat

☐ (480810-E)


cushion to reduce the seat depth (= 15 cm) and the height of the back (= 25 cm)

☐ (480700-E)

You can offer to very small children a position where they are able to put their feet on the seat.



abduction block, removable

☐ (480400-E)


XL-Kit (for a crown height > 62 cm)

☐ (480300-E)

- longer belt and other upholstery
- longer backrest cushion (2 cm thickness)
- if applicable an adapter plate
- instructions

☐ retrofit kit*
(480350-E)

*type serial number please:



further options

protection mat
(to protect original seats of the car)

☐ (480780-E)

belt to connect the chest straps
(to improve the position of the belt in front of the breast)

☐ (480640-E)


security pad for the buckle

☐ with snap fastener
(480650-E)

☐ with magnetic lock
(480655-E)


notes

.....

.....

.....

.....

.....

.....