

Auf der Gröb 2 c  
83064 Raubling  
Germany  
Tel.: +49 (0)8035-9580448  
Fax: +49 (0)8035-9580446  
info@hernik.de • www.hernik.de



invoice / delivery address

Kom.:

**Order**

**Offer**

# Rehab Car Seat **KIDSFLEX R129** NEW! Edition

**Exportprices** **Release: 2024**

## BASIC

	size 1	size 2
altura total	74 cm	80 cm
crown height	67 cm	73 cm
back height	37 - 51 cm	41 - 57 cm
seat width	31 cm	35 cm
seat depth (adjustable)	28 - 33 cm	34 - 39 cm
lateral support height (small)	26 - 36 cm	26 - 36 cm
lateral support height (large)	30 - 40 cm	30 - 40 cm
distance of lateral supports	16 - 28 cm	19 - 32 cm
weight	8,2 kg	8,9 kg
capacity	32 kg	36 kg

**common information (standard):** NEW!

- seat depth adjustable
- headrest included (longer backrest)
- upholstery & cover
- backrest angle adjustable
- 5-point-harness with upholstery
- frame to fix a footrest, clamp bow or ISOFIX
- Crash tested according to **R129/03**
- **ISOFIX adapter** included
- **headrest KICK-BACK** included
- **trunk supports** small, fix, included
- **shoulder support** included



size 1  
Black  
(300000-E3)



size 2  
Black  
(310000-E3)

Technische Änderungen und Druckfehler vorbehalten.

**belt guide\*:**

- LEFT
- RIGHT

\*please choose ONE of them

wider hip cushions (SW -4 cm)  
(300850-E)

**fittings for the different needs**

head support „Kick-Back“  
with aluminum frame  
(free of charge)



head fixation  
(300630-E)



leg support

left  
(300350-E)

right  
(300355-E)



abduction block, removable

size 1  
(300400-E)

size 2  
(310400-E)

size 2XL  
(310400-E-XL)



foot rest, height and angle adjustable

size 1  
(300450-E)

size 2 und 2XL  
(310450-E)


swing up  
(300452-E)



foot roll, height and angle adjustable

size 1  
(300455-E)

size 2 und 2XL  
(310455-E)



trunk supports, small  
(height: 11 cm, depth: 14 cm)

fix  
(free of charge)

swing away  
(300505-E3)



trunk supports, large  
(height: 21 cm, depth: 16 cm)

fix  
(300510-E3)

swing away  
(300515-E3)



shoulder supports  
(free of charge)




**fittings frame**

swivel base

size 1, 2, 2XL  
(300200-E)

RIGHT

LEFT



fastener

(independent of the direction of rotation)

inside  
(standard)

outside  
(300220-E)

angular adjustment (0° - 10°)

size 1  
(300100-E)

size 2  
(310100-E)

size 2XL  
(310100-E-XL)



extra angular adjustment (10° - 20°)

Gr. 1  
(300105-E)

Gr. 2  
(310105-E)

Gr. 2XL  
(310105-E-XL)



Sandwich-Base (mounted under the base board) (300790-E)

Useful in the following situations:

- with very sporty shaped car seats, in order to compensate for the recess there. This results in a more stable position and more freedom of movement when using the swivel adapter.
- with car Isofix brackets that are very far inside to allow nevertheless contact with the car's backrest and the IPAI-NXT (the ISOFIX adapter sits in the sandwich base)
- approx. 1.5 kg additional weight



