

Auf der Gröb 2 c  
83064 Raubling  
Germany  
Tel.: +49 (0)8035-9580448  
Fax: +49 (0)8035-9580446  
info@hernik.de • www.hernik.de

invoice / delivery address

Reference:

**Order**

**Offer**

# Rehab Car Seat HERO XXL

**Exportprices** **Release: 2024**

**BASIC**

**Details:**

HERO-XL	
SEAT WIDTH (pelvis)	24 cm
SEAT WIDTH (knee)	27 cm
SEAT DEPTH	25 - 34 cm
BACK HEIGHT	28 - 45 cm
CROWN HEIGHT	48 - 68 cm
HEIGHT LATERALS	13 - 20 cm
total height	65 - 82 cm
total width	43 cm
total depth	55 cm
capacity	36 kg

**common information (standard):**

- homologated according to ECE-R 129/03
- 5-point-harness with upholstery
- height adjustable head rest
- additional upholstery for the seat, back and head rest
- large seat inclination
- tie-down straps
- side impact protection
- covers washable (30 °C)



**Note:** ISOFIX latches in the car are NOT necessary!!!

Technische Änderungen und Druckfehler vorbehalten.



**accessories for stabilization**

**MOUNTING FRAME**  
(necessary for the options on the right!)  
(490100-E)

**abduction block**  
(490400-E)



**seatdepth adapter**  
(490850-E)




**further accessories**

**security tray**  
(490600-E)



**side supports made of aluminum, fully adjustable**  
(490800-E)



**security pad for the buckle with snap fastener**  
(490650-E)



**belt to connect the chest straps**  
(490640-E)



**protection mat**  
(to protect the seats of the car)  
(490780-E)



**head fixation**  
(in case of poor head control)  
(490630-E)



**notes**

**Tether straps**  
(490710-E)



**Attachment loops, 70 mm**  
(490720-E)



**Lower Tether Anchorages**  
120 mm  
(490730-E)



**Lower Tether Anchorages**  
160 mm  
(490740-E)



Notes section with horizontal dashed lines for writing.