

invoice / delivery address

Reference:

Order

Offer

Rehab Car Seat **IPAI^{MAX}**

Exportprices **Release: 2025**


BASIC

Details:

	IPAI-MAX WITHOUT extra cushion	IPAI-MAX WITH extra cushion
SEAT WIDTH (pelvis)	32 cm	26 cm
SEAT WIDTH (knee)	38 cm	38 cm
SEAT DEPTH	33 cm	30 cm
BACK HEIGHT	41 - 56 cm	38 - 53 cm
CROWN HEIGHT	59 - 76 cm	56 - 73 cm
HEIGHT LATERALS	24 - 36 cm	21 - 33 cm
DISTANCE LATERALS	15 - 28 cm	
angle	80° - 120°	
total height	71 - 88 cm	
total width	44 cm	
total depth	44 cm	
weight	from 8,0 kg	

common information (standard):

- crash tested according to ECE-R129/03
- 5-point-harness with upholstery, height adjustable
- belt mechanics installed into the front of the seat
- solid aluminium board, including the mechanism to recline the backrest stepless
- head rest height adjustable, including belt guides
- additional upholstery (20 mm) for the back & headrest
- ISOFIX adapter or clamp bow integrated in base board



IPAI-MAX red-black
(450000-E1)



IPAI-MAX grey-black
(450000-E2)



Standard belt position for belt buckle in the seat cushion
(19 cm from back cushion)

Additional belt position for belt buckle in the seat cushion
(14 cm from back cushion)



accessory that makes the base board NECESSARY


BASE BOARD

Necessary for the options below;
incl. ISOFIX adapter or clamp bow



ISOFIX

(450000-BEI)



(450000-BEK)



swivel base

(450200-E)

right

left



angular adjustment (5° - 20°)
(450100-E)



Conversion set

from ISOFIX adapter
to clamp bow
(450020-E)

from clamp bow to
ISOFIX adapter
(450030-E)

foot rest

(450450-E)




swing up*
(450452-E)



*450450-E necessary!


outrigger*
(450453-E)



*450450-E necessary!

foot roll

(450455-E)



further options

protection mat

(to protect the original seats of the car)

(450780-E)



Security pad for the buckle

cover with snap fastener
(450650-E)



with HERNIK Security Key*
(450655-E)



*incl. belt cutter



extra cushions

<input type="checkbox"/> back cushion (450700-E)	<input type="checkbox"/> seat cushion (450720-E)	<input type="checkbox"/> head upholstery, left & right (450740-E)	<input type="checkbox"/> head upholstery, behind (450750-E)	<input type="checkbox"/> upper body upholstery, pair (450710-E)	<input type="checkbox"/> hip upholstery, pair (450730-E)
reduces seatdepth ~3 cm	reduces back height ~3 cm	head rest gets more tight ~5 cm	head comes 2 cm forward	reduction of back width ~6 cm	reduction of seat width ~6 cm

fittings for stabilization

<input type="checkbox"/> abduction block, removable (450400-E)	<input type="checkbox"/> head fixation (450630-E)	<input type="checkbox"/> security tray (450600-E)	<input type="checkbox"/> side supports, pair, adjust. in height, depth, width (450800-E)
<input type="checkbox"/> seatdepth adapter (450850-E)	<input type="checkbox"/> lateral support (450810-E)	<input type="checkbox"/> belt to connect the chest straps (450640-E)	

Notes

.....

.....

.....

.....