

invoice / delivery address

Reference:

Order

Offer

Rehab Car Seat **IPAI** **MINI**

Exportprices **Release: 2025**


BASIC

Details:

	IPAI-MINI
SEAT WIDTH (pelvis)	25 cm
SEAT WIDTH (knee)	31 cm
SEAT DEPTH	26 cm
BACK HEIGHT	30 - 43 cm
CROWN HEIGHT	52 - 65 cm
HEIGHT LATERALS	17 - 31 cm
DISTANCE LATERALS	14 - 23 cm
angle	80° - 120°
total height	66 - 89 cm
total width	44 cm
total depth	44 cm
weight	ab 8 kg

common information (standard):

- crash tested according to ECE-R129/03
- 5-point-harness with upholstery, height adjustable
- belt mechanics installed into the front of the seat
- solid aluminium board, including the mechanism to recline the backrest stepless
- head rest height adjustable, including belt guides
- additional upholstery (20 mm) for the back & headrest
- ISOFIX adapter or clamp bow integrated in base board



IPAI-MINI red-black
(480000-E1)



IPAI-MINI grey-black
(480000-E2)




accessory that makes the base board NECESSARY

BASE BOARD

Necessary for the options below!



(480000-BEI)



(450000-BEK)



swivel base

(480200-E)

right

left



angular adjustment
(5° - 20°)
(480100-E)



Conversion set

from ISOFIX adapter
to clampbow
(480020-E)

from clampbow
to ISOFIX adapter
(480030-E)

foot rest

(480450-E)




swing up*
(480452-E)



*480450-E necessary!


outrigger*
(480453-E)



*480450-E necessary!

foot roll

(480455-E)



further options

protection mat

(to protect the original seats of the car)

(480780-E)



security pad for the buckle

with snap fastener
(480650-E)




with magnetic lock
(480655-E)




fittings for stabilization

abduction block, removable
(480400-E)
[€ 275,00]

A photograph of a grey, rectangular abduction block being attached to the side of a car seat's seat pan.


head fixation
(480630-E)
[€ 375,00]

A photograph of a black head fixation device with a white strap, designed to secure a patient's head.

security tray
(480600-E)
[€ 465,00]

A photograph of a black, rectangular security tray mounted on the front of a car seat.


side supports, pair, adjust.
in height, deep, width
(480800-E), [€ 425,00]

A diagram showing a pair of black side supports on a car seat. Red arrows indicate that the supports can be adjusted in height, depth, and width.

seatdepth adapter
(480850-E)
[€ 315,00]

A photograph of a seatdepth adapter. A red arrow points to the adapter's depth, with the text "+ 2 - 12 cm" next to it.

lateral support
(480810-E)
[€ 325,00]

A photograph of a car seat with two black lateral support pads positioned on either side of the seat pan.

belt to connect the chest
straps
(480640-E), [€ 125,00]

A photograph of a black webbing belt with a plastic buckle, used to connect chest straps.

Notes

A large rectangular area containing horizontal dotted lines for taking notes.